

# Color Choice

- Automatic Color measurement
- Eliminates operator bias
- Performs the most widely used Color scales in the Petroleum industry in seconds
- Color tests performed ASTM D1500, Saybolt D156, D6045, D1209, D5386, ISO 2049, IP 196
- State of the Art Spectrophotometer
- No moving parts or filters to degrade
- On site calibration
  - No internet connection needed
- UL, TUV, CE and CSA Certified
- Simple touch screen
  - Option for keyboard and mouse
- Store and print hundreds of results
- Easy to use, calibrate and maintain
- Output to LIMS, printer, PC or external storage device
- Proudly made in the USA



Available exclusively from:



**Choice Analytical, Inc.**

*Support You Can Trust*

For product and sales information contact: [sales@choiceanalytical.com](mailto:sales@choiceanalytical.com) or 832-271-7448  
[www.choiceanalytical.com](http://www.choiceanalytical.com)